

Fig. 1

Figure 2

SEQ ID NO	Mutation	Exemplary Pool	3' blocking group	Oligo type (Arm)	Sequence (5'-3')	ϵ_{calc} M ² cm ⁻¹
SEQ ID NO: 296	2789+5G>A	1	none	invader	TTTGGTTGCTGTGGCTCTCTGGAAAGTGAT	330800
SEQ ID NO: 297	2789+5G>A	1	hex	probe/DM	CGCCCGGAGGATATCCATGTCTCTATTGTG	306500
SEQ ID NO: 298	2789+5G>A	1	none	synthetic target	CAATCTCACATAAGACATGGAAATATTCATCTCCAAAGGCGCACGACGACACAAA	667000
SEQ ID NO: 299	R1162X	1	none	invader	GTTTACCTCTGTGGCATGTGATGAATCTAAAGACTCT	428000
SEQ ID NO: 300	R1162X	1	hex	probe/DM	CGCCCGGAGGAGCTCACAGTCGC	253000
SEQ ID NO: 301	R1162X	1	none	synthetic target	TCAGATGCCATCTGTGAGCTGAGCTTTAAGTTCATTCAGATGCCACAGAAAGGTTAAAC	659000
SEQ ID NO: 302	R347P	4	none	invader	CAGGAAATTTGCGAGTGCACCCCATGT	306600
SEQ ID NO: 303	R347P	4	hex	probe/ER24	ACGACGCGGAGGCGCAAGAAATGCGAG	318200
SEQ ID NO: 304	R347P	4	none	synthetic target	CTCATCTGCATGTCTGCCATCGCGGTCACTCGGCAATTTCCCTGGG	488000
SEQ ID NO: 305	1898+1G>A	1	none	invader	GACTCTCTTTGGATACCTAGATGTTTAAACAGAAAAGAAATATTTGAAAGT	619800
SEQ ID NO: 306	1898+1G>A	1	hex	probe/DM	CGCCCGGAGGATGTCTTTGAATACCTTACTAT	386500
SEQ ID NO: 307	1898+1G>A	1	none	synthetic target	ATAAGTAAGTATCAAGACATCACTTCAAAATATTTCTTTCTGTATAAACATCTAGGTATCCAAAAGGAGAGTTC	904800
SEQ ID NO: 308	2184delA	4	none	invader	CCCCAACTCTCCAGTCGTTTAAAGATATTTTTC	393000
SEQ ID NO: 309	2184delA	4	hex	probe/DM	CGCCCGGAGGTTCTGTCCAGGAGACA	305200
SEQ ID NO: 310	del507	1	none	invader	GCCTTGATGACGCTCTGTATCTATTCATCATCATAGGAACACCAAT	509300
SEQ ID NO: 311	del507	1	hex	probe/DM	CGCCCGGAGGATATTTCTTTAATGTGGCC	345200
SEQ ID NO: 312	del507	1	none	synthetic target	GCCTTGGCCACATTAAGAAATATCATTGGTGTGTCTGTATGATGAATAGATACAGAAGCGTCATCAAAAGCATCGC	866800
SEQ ID NO: 313	G85E	4	none	invader	GCCTCTGGCGGATTTCTTGGAGATTAAGTCTATGT	409100
SEQ ID NO: 314	G85E	4	hex	probe/ER24	ACGACGCGGAGGAATCTTTTATATTAGGGTAAAG	431700
SEQ ID NO: 315	G85E	4	none	synthetic target	AGATCTATGCCCTTAATATAAAGATTTTCATAGACATAAATCTCCAGAAAAAACAATCGCCGGAAGGGGCAITTA	869300
SEQ ID NO: 316	R117H	3	none	invader	AATCATGCTCTCTATGATCGCCGGATAGCAAGAGGAGAACT	443800
SEQ ID NO: 317	R117H	3	hex	probe/DM	CGCCCGGAGGACTCTATCGCAATTATCT	304200
SEQ ID NO: 318	R117H	3	none	synthetic target	ATCGCTAGATAAATCGCGATAGAGTGTCTCTCTGTATCGGGGTATAGGAAGCTATGATTT	681700
SEQ ID NO: 319	R560T	1	none	invader	CATGAATGACATTTACAGCAATGCTCTAGACCAATATAGTATTCACT	595000
SEQ ID NO: 320	R560T	1	hex	probe/ER24	ACGACGCGGAGGTCTGTAAGAAATTTGCT	378100
SEQ ID NO: 321	R560T	1	none	synthetic target	CAACGACGAGAAATTTCTTGAACAGTGTGAATACTAATATTGGCTAGCAAGCAATTTGCTGTAAATGTCAATGTATAAA	945400
SEQ ID NO: 322	3120+1G>A	2	none	invader	GCAATTTGGATGACCTCTGCGCTCTTACCATTATTGACTTTCATCCAGT	486000
SEQ ID NO: 323	3120+1G>A	2	hex	probe/DM	CGCCCGGAGGATGTATAAAATAGTACCGTTAA	397500
SEQ ID NO: 324	3120+1G>A	2	none	synthetic target	AGACATCTAATCCGATCTATTTTATCATCTCTGATGAAGTCAAATATGTAAGGCGAGAGGTCATCCAAAATGCTATATC	984000
SEQ ID NO: 325	3659delC	2	none	invader	GAGATGTGCCATCTGTGAAGTTTGGTGTGACTTT	372300
SEQ ID NO: 326	3659delC	2	hex	probe/DM	CGCCCGGAGGTAGTTACCTCTCTGTGG	302800
SEQ ID NO: 327	3659delC	2	none	synthetic target	CATGCCACAGAAAGTTAAACTCAAGTCAAGTACCAACCAACCATCAAGAAATGGCCCAACTCTC	679800
SEQ ID NO: 328	A455E	1	none	invader	CTCTGAAGATATTAAATCTCAAGTAGAAAGGACGAGGAGCTGTGGT	531000
SEQ ID NO: 329	A455E	1	hex	probe/ER24	ACGACGCGGAGAGTTGCTGGATCCA	298100
SEQ ID NO: 330	A455E	1	none	synthetic target	CCAGTGGATCCACCAACTCCACCAACTGCTCTCTTCTATCTTGAATTAATATCTTTCCAGG	661000
SEQ ID NO: 331	1078delT	2	none	invader	AGTGCATCGGAAGCAGATATAAACACCAAT	413500
SEQ ID NO: 332	1078delT	2	hex	probe/DM	CGCCCGGAGGAGCCCTGAGAAGAGAA	355400
SEQ ID NO: 333	1078delT	2	none	synthetic target	AGCCTCTCTCTCAGGGTCTGTGGTGTGTTTATCTGTCTCTCCCTATGCACT	533300
SEQ ID NO: 334	G551D	2	none	invader	GCAGAGAAAGACAATATGATCTTGTGGAAGGTGGAAACACATGAGTGGAGT	628200
SEQ ID NO: 335	G551D	2	hex	probe/DM	CGCCCGGAGGATCAACGACGACGAATTTCT	343800
SEQ ID NO: 336	G551D	2	none	synthetic target	CTCTGTAAGAAATTTCTGCTGTGTTCTCTCCACTCAGTGTGATTCACCTCTCTCCAAGCACTATATTGCTCTCTGCAAACTT	883100
SEQ ID NO: 337	I148T	1	none	invader	AAATCAACATAACTAGCTATCTCATCTGCATTCCAT	432400
SEQ ID NO: 338	I148T	1	hex	probe/ER24	ACGACGCGGAGGTGTGATGAAGGCGCAA	350200
SEQ ID NO: 339	I148T	1	none	synthetic target	CCATTTTGGCTCTCATCACTGGAATGCAGATGAGTAAGTACGTATGTTAGTTGATTT	643100
SEQ ID NO: 340	N1303K	2	none	invader	CCATATTTCTTGATCACTGCTGCTTATGGGATGCCAAT	414700
SEQ ID NO: 341	N1303K	2	hex	probe/DM	CGCCCGGAGGCTTTTCTGTAATTTGTCAGAAAAT	391200
SEQ ID NO: 342	N1303K	2	none	synthetic target	ATTATTTTCTCGAACATTTAGAAAAGTTGGATCCCTTATGAACAGTGGAGTGATCAAGAAATATGGAAG	867100
SEQ ID NO: 343	711+1G>T	2	none	invader	GCCTTCCAGTTGTATATTTATAACAAATAGTGGCTTAAGAAATTAAGAAATTAACAATAGTACATT	
SEQ ID NO: 344	711+1G>T	2	hex	probe/DM	CGCCCGGAGGAATTCATAAATTTCTCAGGT	
SEQ ID NO: 345	711+1G>T	2	none	synthetic target	ACCTGAAACAAATTCATGAATATGTACACTATGATTTAATCTTTTAGGCACTATGTATAAATATATCAACTTGGAAAGGC	927000
SEQ ID NO: 346	1717-1G>A	3	none	invader	GCCCTTCAAATTCAGATTCAGGATCACTAATAAAGTGACTCTTAATTTCTTTTGGTAAT	685000
SEQ ID NO: 347	1717-1G>A	3	hex	probe/DM	CGCTTTCAGGAGTGTATGACATCTCCAAAGTTGC	294500
SEQ ID NO: 348	1717-1G>A	3	none	synthetic target	CTCTGCAAACTTGGAGTACTCTCTGCTGTGTTCTCTCCACTCAGTGTGATTCACCTCTTATAGTATGCTCAATCTGAAATTTGAAAGGCACATC	1010000
SEQ ID NO: 349	W1282X	3	none	invader	GCTCACTGTGTATGCTCTCCAAAGTTCTCTA	345000
SEQ ID NO: 350	W1282X	3	hex	probe/DM	CGCCCGGAGGTACTGTTCGCAAGTTATG	327800
SEQ ID NO: 351	W1282X	3	none	synthetic target	GATTCAATACTTCTGCAACAGTGAAGGAAGCCCTTGGAGTGTATACCAAGGTGAGCAA	683000
SEQ ID NO: 352	3849+10KbC-T	2	none	invader	CAAGAGTCTCCATCTGTTCGAGTATTAATAATGGA	390000
SEQ ID NO: 353	3849+10KbC-T	2	hex	probe/DM	CGCCCGGAGGAGTGAAGACACCCCTGAAA	327400
SEQ ID NO: 354	3849+10KbC-T	2	none	synthetic target	TTCCTTCAGGCGTGTATGATCTACCACTTTAATACTGCAACAGATGAAGACTCTTGT	
SEQ ID NO: 355	R553X	4	hex	probe/DM	CATTACAGCAAAATGCTCTGACAGCAATATAGTATTATCCACTTTCGTAAGAAATTTCTTGCTG	
SEQ ID NO: 356	R553X	4	none	invader	CGCCCGGAGGCTGACTGCTCCACTCAT	

FIGURE 3

SEQ ID NO	Mutation	Oligo Name	Oligo type	Sequence (5'-3')	ε260 M-1 cm-1	Length
SEQ ID NO:1	A455E	1297-83-03	Invader	CCTGAAGATATTAATTTCAAGATAGAAAGGACAGTTGTTGGT	465900	45
SEQ ID NO:2	A455E	1657-96-01	WT ER38	AGGCCACGGACGGGTTGCTGGATC-hex	235400	29
SEQ ID NO:3	A455E	1297-31-05	MT ER24	ACGGACGGAGAGTGTGGATCCCA-hex	266300	31
SEQ ID NO:4	A455E	1618-54-14	WT Syn Target (ST)	CCAGTGGATCCAGCAACGCCCAACACTGCTCTTCTATCTTGAATTAATATCTTTCAGGT	600200	64
SEQ ID NO:5	A455E	1618-54-18	Mut Syn Target (ST)	CCAGTGGATCCAGCAACCTCCCAACACTGCTCTTCTATCTTGAATTAATATCTTTCAGGT	598800	64
SEQ ID NO:6	A455E	23-204	FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:7	A455E	23-755	Red FRET/FAM Stem-MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:8	3659delC	1618-60-10	Invader	AGTTGATTTGACATGCCAACCAAGGTAAACCTAT	343500	34
SEQ ID NO:9	3659delC	1618-60-12	WT ER24	ACGGACGGAGCGCAAGTCAACCAACCA-hex	290200	33
SEQ ID NO:10	3659delC	1618-60-17	MT ER38	AGGCCACGGACCGCAAGTCAACCAACCACTAC-hex	307600	35
SEQ ID NO:11	3659delC	1618-60-14	WT ST	TCCTGTATGGTTGGTGGTACTTGGTGGTTCCTCTCTGTGGCATGTCAATGAAC	531400	57
SEQ ID NO:12	3659delC	1618-60-18	MT ST	TCCTGTATGGTTGGTGGTACTTGGTGGTTCCTCTCTGTGGCATGTCAATGAAC	521300	56
SEQ ID NO:13	3659delC	23-210	FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:14	3659delC	23-211	Red FRET - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:15	G85E	1614-57-10	Invader	GCCCTTCGGGATGTTTTTCTGGAGATTATGTTCTATGT	374700	41
SEQ ID NO:16	G85E	1614-57-13	WT ER24	ACGGACGGGAGGAATCTTTTATATTAGGGTAAAG-hex	372800	41
SEQ ID NO:17	G85E	1614-57-18	MT ER38	AGGCCACGGACGAATCTTTTATATTAGGGTAAAG-hex	381500	42
SEQ ID NO:18	G85E	1614-57-12	WT ST	AGATCTTACCCTTAATAAAGATTCATAGAACATAAATCTCCAGAAAAACATCGCGAAGGGCAAT	742300	73
SEQ ID NO:19	G85E	1614-57-16	MT ST	AGATCTTACCCTTAATAAAGATTCATAGAACATAAATCTCCAGAAAAACATCGCGAAGGGCAAT	743200	73
SEQ ID NO:20	G85E	1055-48-09	FAM FRET-WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:21	G85E	1230-33-01	Red FRET/FAM Stem-MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:22	N1303K	16145404	WT Probe ER24	ACGGACGGGAGGTTTTTCTAAATGTTCCAGAAAAA-hex	369300	41
SEQ ID NO:23	N1303K	16145406	MT Probe DM	CGCGCGGAGGTTTTTCTAAATGTTCCAGAAAAA-hex	334500	39
SEQ ID NO:24	N1303K	16145401	Invader	TCTTGATCACTCCACTGTTCTAGGGATCCAA	308600	34
SEQ ID NO:25	N1303K	16145403	WT ST	TATTTTCTGGAACATTTAGAAAAAAGTGGATCCCTATGAACAGTGAGTGATCAAGAAAT	650400	65
SEQ ID NO:26	N1303K	16145407	MT ST	TATTTTCTGGAACATTTAGAAAAAAGTGGATCCCTATGAACAGTGAGTGATCAAGAAAT	654900	65
SEQ ID NO:27	N1303K	23-428	DM FAM FRET MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-ggc-hex		
SEQ ID NO:28	N1303K	23-384	ER24 RED FRET WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-ct-hex	451000	46
SEQ ID NO:29	D1270N	1645-18-10	Invader	TACGTCTTTTGGAGACTACTGAACACTGGATTCCTCTCAGTGTTCAGTGTCTCAAAAAAGCTGATAT	257300	31
SEQ ID NO:30	D1270N	1665-93-01	WT DM	CGCGCGGAGGATGTTGTTCTGGGA-hex	275600	32
SEQ ID NO:31	D1270N	1645-18-17	MT Probe (ER24)	ACGGACGGGAGGATGTTGTTCTGGGA-hex	659200	68
SEQ ID NO:32	D1270N	1645-18-12	WT ST	GAATCCCAAGACACACCATGATCTGGATTCCTCTCAGTGTTCAGTGTCTCAAAAAAGCTGATAT	659500	68
SEQ ID NO:33	D1270N	1645-18-16	MT ST			
SEQ ID NO:34	D1270N	23-746	DM Red FRET/FAM stem- WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-ggc-hex		
SEQ ID NO:35	D1270N	23-210	ER24 FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:36	R560T	1614-55-01	Invader	GACATTTACAGCAATGCTTCTGACCACTAATATTGTTTCACT	452500	46
SEQ ID NO:37	R560T	1614-55-04	WT probe (ER24)	ACGGACGGGAGGTTGCTAAAGAAATCTTGTCT-hex	308700	37
SEQ ID NO:38	R560T	1614-55-06	MT probe (DM)	CGCGCGGAGGTTGCTAAAGAAATCTTGTCT-hex	283300	35
SEQ ID NO:39	R560T	1614-55-03	WT ST	AACGACGAAGATTTCTTTAGCAAGCTGAATACTAATATTGGTCTAGCAAGCATTTGCTGTAATGTCT	706600	71
SEQ ID NO:40	R560T	1614-55-07	MT ST	AACGACGAAGATTTCTTTAGCAAGCTGAATACTAATATTGGTCTAGCAAGCATTTGCTGTAATGTCT	703300	71
SEQ ID NO:41	R560T	23-210	ER24 FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:42	R560T	23-205	DM Red FRET - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-ggc-hex		
SEQ ID NO:43	621+1G>T	1618-48-10	Invader	CCTCATCACATTTGAAATGAGATGAGATAGCTATGTTAGTTGATTTAATAAGC	588400	59
SEQ ID NO:44	621+1G>T	1665-71-02	WT ER24	ACGGACGGGAGGTAATCTCTCTGCAC-hex	277800	33
SEQ ID NO:45	621+1G>T	1618-48-15	MT Probe (DM)	CGCGCGGAGGTAATCTCTCTGCACAGG-hex	267700	35
SEQ ID NO:46	621+1G>T	1618-48-14	WT ST	GGGGCTGTGAAGGAATTAACCTCTTATAAATAAACTAAACATAGCTATTCATCTGCATTCCTCAATGTGAAGG	825700	85
SEQ ID NO:47	621+1G>T	1618-48-18	MT ST	GGGGCTGTGAAGGAATTAACCTCTTATAAATAAACTAAACATAGCTATTCATCTGCATTCCTCAATGTGAAGG	830500	85
SEQ ID NO:48	621+1G>T	23-755	ER24 Red FRET w/FAM stem- WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:49	621+1G>T	23-428	DM FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-ggc-hex		
SEQ ID NO:50	1717-1G>A	1548-76-07	Invader	GATGTCCTTTCAATTCAGATGAGCATCTAAAGTGACTCTCTAATTTCTATTTTGGTAACTC	644800	67
SEQ ID NO:51	1717-1G>A	1657-96-02	WT probe (DM)	CGCGCGGAGGAGACATCTCCAAGTTG-hex	286500	32
SEQ ID NO:52	1717-1G>A	1657-40-02	MT probe (ER24)	ACGGACGGGAGGAGACATCTCCAAGTTG-hex	298500	34
SEQ ID NO:53	1717-1G>A	1618-56-14	WT ST	CTCTGCAAACTGGAGATGCTATTACCAAAATAGAGATGACCTTTAGTATGCTCAATCTGAAATTTGAAAGG	863600	87
SEQ ID NO:54	1717-1G>A	1618-56-18	MT ST	CTCTGCAAACTGGAGATGCTATTACCAAAATAGAGATGACCTTTAGTATGCTCAATCTGAAATTTGAAAGG	864500	87

FIGURE 3

SEQ ID NO	Mutation	Oligo Name	Oligo type	Sequence (5'-3')	ε260 M-1 cm-1	Length
SEQ ID NO:55	1717-TG>A	23-746	DM Red FRET w/FAM stem - WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-gg-hex		
SEQ ID NO:56	1717-TG>A	23-210	ER24 FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:57	1078delT	1614-80-01	Invader	CCCTTGATTAGTCATAGGGAACACAGATAAAACACACAT	434600	43
SEQ ID NO:58	1078delT	1614-80-02	WT Probe DM	CGCGCGAGGAAGAACCCCTGAGAGAAGA-hex	298900	29
SEQ ID NO:59	1078delT	1614-80-03	MT Probe ER24	ACGACGCGGAGAGACCCCTGAGAGAAGA-hex	329000	31
SEQ ID NO:60	1078delT	1614-80-08	WT ST	AGCCTTCTCTCTCAGGGTCTTGTTGTTTATCTGTGCTTCCCTATGCACTAATCAAGGAA	623300	68
SEQ ID NO:61	1078delT	1614-80-07	MT ST	AGCCTTCTCTCTCAGGGTCTTGTTGTTTATCTGTGCTTCCCTATGCACTAATCAAGGAA	615200	67
SEQ ID NO:62	1078delT	23-428	DM FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:63	1078delT	23-755	ER24 Red FRET w/FAM Stem - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:64	R347P	1618-62-01	Invader	CAGGGAATTCGCCAGTAGCCGCAATGT	271100	28
SEQ ID NO:65	R347P	1697-40-04	WT Probe (DM)	CGCGCGGCGGCGCAACAATGCAGAAT-hex	275800	32
SEQ ID NO:66	R347P	1696-93-05	MT ER24	ACGACGCGGAGGCGCAACAATGCAGAAT-hex	266700	30
SEQ ID NO:67	R347P	1618-62-05	WT ST	CTCATCTGCATGTTCTGCCATGGCGGTCACTGGCAATTCCTCGGT	452500	51
SEQ ID NO:68	R347P	1618-62-09	MT ST	CTCATCTGCATGTTCTGCCATGGCGGTCACTGGCAATTCCTCGGT	450200	51
SEQ ID NO:69	R347P	23-394	ER24 Red FRET - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:70	R347P	23-428	DM FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:71	2184delA	1532-85-01	WT probe (DM)	CGCGCGAGGGTCTTCTGTCACGAGAG-hex	257300	32
SEQ ID NO:72	2184delA	1582-25-02	Invader	CGCGCGGCGGAGTCTTCTGTCACGAGAG-hex	295500	35
SEQ ID NO:73	2184delA	1532-85-04	WT probe (DM)	CGCGCGGCGGAGTCTTCTGTCACGAGAG-hex	418200	46
SEQ ID NO:74	2184delA	1680-88-17	plasmid Het target	CTTCTTTTCCCCAACTCTCCAGTCTGTTTAAAGATTGTTTA		
SEQ ID NO:75	2184delA	23-210	ER24 FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:76	2184delA	23-746	DM Red FRET w/FAM stem - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:77	V520F	16186602	WT Probe (DM)	CGCGCGGCGGCTCTGTATCTATATCATCAT-hex	304500	33
SEQ ID NO:78	V520F	16186607	MT Probe (ER24)	ACGAGCGGCGGAGAGTCTGTATCTATATCATCAT-hex	348600	36
SEQ ID NO:79	V520F	16186607	Invader	CACATAGTTTCTTACTCTCTCTAGTTGGCATGCTTTGATGAT	386600	42
SEQ ID NO:80	V520F	16186605	WT Synthetic Target	TCCTATGATGATATAGATACAGAGCTTCATCAAGCATGCCAACTAGAGAGGTAAGAACTATGTGAAT	747300	72
SEQ ID NO:81	V520F	16186609	MT Synthetic Target	TCCTATGATGATATAGATACAGAGCTTCATCAAGCATGCCAACTAGAGAGGTAAGAACTATGTGAAT	744100	72
SEQ ID NO:82	V520F	23-428	DM FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:83	V520F	23-394	ER24 Red FRET - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:84	R347H	1618-68-10	Invader	GGAAATATTCACCATCTCATCTGCTGCTGCTGCT	370700	40
SEQ ID NO:85	R347H	1618-68-12	WT Probe (ER24)	ACGAGCGGCGGAGCATGGCGGTAC-hex	242300	29
SEQ ID NO:86	R347H	1665-93-04	MT ER38	AGCCACGACGAGAGCAATCTTTTAAACAGACTG-hex	242300	29
SEQ ID NO:87	R347H	1618-68-14	WT ST	GCCGAGTGACCGCATGCGCAGAACATGCAGAATGAGATGGTGGTGAATTTTCCT	643800	58
SEQ ID NO:88	R347H	1618-68-18	MT ST	GCCGAGTGACCGCATGCGCAGAACATGCAGAATGAGATGGTGGTGAATTTTCCT	645100	58
SEQ ID NO:89	R347H	23-210	ER24 FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:90	R347H	23-752	ER38 Red FRET w/FAM stem - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:91	2183AA>G	1645-26-10	Invader AS	GAGATGCTCTCTGCTCTGAGACAGCAACAAAT	335200	34
SEQ ID NO:92	2183AA>G	1667-19-08	WT ER24	ACGAGCGGAGAGCAATCTTTTAAACAGACTG-hex	331900	37
SEQ ID NO:93	2183AA>G	1667-19-10	MT DM	CGCGCGGAGGCAATCTTTTAAACAGACTG-hex	286500	34
SEQ ID NO:94	2183AA>G	1645-26-12	WT ST	CCCCAACTCTCCAGTCTGTTTAAAGATTGTTTGTCTGTCGAGACAGGACATCTCCTTCT	654900	71
SEQ ID NO:95	2183AA>G	1645-26-16	MT ST	CCCCAACTCTCCAGTCTGTTTAAAGATTGTTTGTCTGTCGAGACAGGACATCTCCTTCT	644100	70
SEQ ID NO:96	2183AA>G	23-755	ER24 Red FRET w/FAM stem - WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:97	2183AA>G	23-428	DM FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:98	R334W	1267-83-04	Invader	CGCAGAACATGCAGAATGAGATGGTGGTGAATTTTCCT	409900	41
SEQ ID NO:99	R334W	1618-64-02	WT DM	CGCGCGGAGGAGGATGATCTCTTGTAT-hex	267600	34
SEQ ID NO:100	R334W	1618-64-07	MT ER24	ACGAGCGGAGAGGATGATCTCTTGTAT-hex	316900	37
SEQ ID NO:101	R334W	1618-64-05	WT ST	GCATTAATCAAGGAATCATCTCTCGGAAATATTCACCACCATCTCATCTGCGATTGTCGCT	624200	66
SEQ ID NO:102	R334W	1618-64-09	MT ST	GCATTAATCAAGGAATCATCTCTCGGAAATATTCACCACCATCTCATCTGCGATTGTCGCT	624500	66
SEQ ID NO:103	R334W	23-746	DM Red FRET w/FAM WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:104	R334W	23-210	ER24 FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gg-hex		
SEQ ID NO:105	R117H	1614-48-10	Invader	AGAATCATAGCTTCTATGACCGGATACCAAGAGGAACT	414400	41
SEQ ID NO:106	R117H	1614-48-11	WT DM	CGCGCGGAGGCTCTATGCGGATTATCTA-hex	277600	34
SEQ ID NO:107	R117H	1614-48-18	MT ER38	AGGCCAGGAGCACTCTATCGGATTATCTAG-hex	317300	37
SEQ ID NO:108	R117H	1614-48-12	WT ST	ATGCCAGTAATAATCGGATAGAGCGTCTCTTGTATTCGGGTCTAGGAAGTATGATTTCTTCT	651200	68

FIGURE 3

SEQ ID NO	Mutation	Oligo Name	Oligo type	Sequence (5'-3')	ε260 M-1 cm-1	Length
SEQ ID NO:109	R117H	1614-48-16	MT ST	ATGCCATAGATAAATCGGATAGAGTGTCTCTCTGTTATCCGGGTATAGGAAGCTATGATCTTCT	653300	68
SEQ ID NO:110	R117H	23-428	DM FAM FRET WT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-hex		
SEQ ID NO:111	R117H	23-211	ER38 Red FRET MT	Z35-tct-X-tcg-gcc-ttt-tcc-ggc-aga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:112	2789+5G>A	1297-83-06	Invader	TTTGGTGTGCTGTGGCTCTCTTGGAAAGTAT	296200	32
SEQ ID NO:113	2789+5G>A	1618-46-13	WT ER38	AGGCCACGACGAGTATCCATGCTCTTATTTGTG-hex	292600	37
SEQ ID NO:114	2789+5G>A	1618-46-16	MT ER24	ACGGACGCGGAGATATCCATGCTCTTATTTGTG-hex	298400	37
SEQ ID NO:115	2789+5G>A	1618-46-14	WT ST	TCTACACAATAGGACATGGAATGACTACTCTTCCAAAGGAGGACAGCACAAACCAAT	554800	56
SEQ ID NO:116	2789+5G>A	1618-46-18	MT ST	TCTACACAATAGGACATGGAATGACTACTCTTCCAAAGGAGGACAGCACAAACCAAT	556700	56
SEQ ID NO:117	2789+5G>A	23-211	ER38 Red FRET WT	Z35-tct-X-tcg-gcc-ttt-tcc-ggc-aga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:118	2789+5G>A	23-210	ER24 FAM FRET MT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex	400100	43
SEQ ID NO:119	394delTT	1665-32-10	Invader	CTTCGGCGGATTTTCTGGAGATTTATGTTCTTATGGAAT	371400	42
SEQ ID NO:120	394delTT	1714-28-02	WT (ER38)	AGGCCACGACGCTTTTATATTTAGGGTAAGGATCT-hex	368000	41
SEQ ID NO:121	394delTT	1680-16-07	MT (ER24)	AGGGACGCGGAGCTTATATTTAGGGTAAGGATCTC-hex	788000	77
SEQ ID NO:122	394delTT	1614-65-12	WT ST	ACAAATGAGATCCTTACCCCTTAAATATAAAGATTTCCATAGAACAATAAATCTCCAGAAAAACATGCGCCGAAGGGC	764000	75
SEQ ID NO:123	394delTT	1614-65-16	MT ST	ACAAATGAGATCCTTACCCCTTAAATATAAAGATTTCCATAGAACAATAAATCTCCAGAAAAACATGCGCCGAAGGGC		
SEQ ID NO:124	394delTT	23-211	ER38 Red FRET - WT	Z35-tct-X-tcg-gcc-ttt-tcc-ggc-aga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:125	394delTT	23-210	ER24 FAM FRET - MT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:126	3849+10KbC>T	1618-44-01	Invader	TGCATGTACCATGAATAGACATTTCTTTCAGGGTGTCTTACTCT	492000	47
SEQ ID NO:127	3849+10KbC>T	1757-50-02	WT DM	CGCGCGAGGGCCATTTTAACTACTGCAACAGA-hex	308400	36
SEQ ID NO:128	3849+10KbC>T	1618-44-07	MT probe (ER24)	ACGGACGCGGAGACCACTTTTAACTACTGCAACAG-hex	372800	37
SEQ ID NO:129	3849+10KbC>T	1618-44-05	WT ST	CCATCTGTTGCAGTATTAATGCGGAGTAAAGACACCCCTGAAAGGAAATGTTCTATTCATGGTACAAATGCAT	715200	72
SEQ ID NO:130	3849+10KbC>T	1618-44-09	MT ST	CCATCTGTTGCAGTATTAATGCGGAGTAAAGACACCCCTGAAAGGAAATGTTCTATTCATGGTACAAATGCAT	717300	72
SEQ ID NO:131	3849+10KbC>T	23-205	DM Red FRET - WT	Z35-tct-X-tcg-gcc-ttt-tcc-ggc-aga-gac-gtc-ggc-ggc-hex		
SEQ ID NO:132	3849+10KbC>T	23-210	ER24 FAM FRET - MT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:133	W1282X	1618-42-01	Invader	GCTCACTGTGGTATCACTTCCAAAGCTTTCTCTA	311700	34
SEQ ID NO:134	W1282X	1618-42-03	WT ER24	ACGGACGCGGAGCACTGTTGCAAGTATTTG-hex	300300	35
SEQ ID NO:135	W1282X	1618-42-06	MT DM	CGCGCGAGGTCACCTGTTGCAAGTATTTG-hex	283200	34
SEQ ID NO:136	W1282X	1618-42-05	WT ST	GATTCATTAACCTTGCAACAGTGGAGGAAAGCCCTTTGGAGTGATACCACAGGTGAGCAAT	602000	60
SEQ ID NO:137	W1282X	1618-42-09	MT ST	GATTCATTAACCTTGCAACAGTGGAGGAAAGCCCTTTGGAGTGATACCACAGGTGAGCAAT	603900	60
SEQ ID NO:138	W1282X	23-210	ER 24 FAM FRET WT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:139	W1282X	23-746	DM Red FRET w/FAM stem MT	Z35-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:140	G542X	1614-49-10	Invader	AATAGGACATCTCCAACTTTGCGAGAAAGACAAATATAGTTCTTC	457000	45
SEQ ID NO:141	G542X	1614-49-14	WT probe (ER38)	AGGCCACGAGCGGAGGAAAGGTGGAAATCAC	296800	29
SEQ ID NO:142	G542X	1614-49-15	MT probe (DM)	CGCGCGAGGTGAGAGGTGGAAATCAC	282400	28
SEQ ID NO:143	G542X	1614-49-12	WT ST	TCAGTGTGATTCACCTTCTCCAAAGAACTATATGCTCTTCTGCAAACTTGAGAGATGCTCTATTT	622000	67
SEQ ID NO:144	G542X	1614-49-16	MT ST	TCAGTGTGATTCACCTTCTCCAAAGAACTATATGCTCTTCTGCAAACTTGAGAGATGCTCTATTT	626800	67
SEQ ID NO:145	G542X	23-752	ER38 Red FRET w/FAM Stem - WT	Z35-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:146	G542X	23-428	DM-FAM FRET - MT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-hex		
SEQ ID NO:147	3120+1G>A	1618-58-10	Invader	GCAATTTTGGATGACCTCTCGCTCTTACCATAITTGACTTTCATCCAGT	496000	49
SEQ ID NO:148	3120+1G>A	1745-68-03	WT ER38	AGGCCACGAGCGGTATGAAAAATAGTACCGTT-hex	348100	38
SEQ ID NO:149	3120+1G>A	1618-58-15	MT probe (DM)	CGCGCGAGGATATGAAAAATAGTACCGTTAA-hex	350500	38
SEQ ID NO:150	3120+1G>A	1618-58-14	WT ST	ATACTAACGGTACTATTTTACATCTGATGAAAGTCAAAATATGTTAGAGGCAAGGTCATCCAAAATTTGCTAT	796400	79
SEQ ID NO:151	3120+1G>A	1618-58-18	MT ST	ATACTAACGGTACTATTTTACATCTGATGAAAGTCAAAATATGTTAGAGGCAAGGTCATCCAAAATTTGCTAT	798300	79
SEQ ID NO:152	3120+1G>A	23-752	ER38 Red w/FAM stem - WT	Z35-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:153	3120+1G>A	23-428	DM FAM FRET - MT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-hex		
SEQ ID NO:154	1148T	1614-50-10	Invader	CTACACCCAGCCATTTTGGCCTTCATCACAA	298400	33
SEQ ID NO:155	1148T	1614-50-11	WT probe (DM)	CGCGCGAGGTGGAATGCAGATGAGAAAT-hex	304400	30
SEQ ID NO:156	1148T	1614-50-17	MT probe (ER24)	ACGGACGCGGAGCTGGAATGCAGATGAGAA-hex	309700	30
SEQ ID NO:157	1148T	1614-50-12	WT ST	TAGCTATTCATCTGCATTCATTCATGTAAGGCAAAATGGCTGGGTGAGGAGT	576700	59
SEQ ID NO:158	1148T	1614-50-16	MT ST	TAGCTATTCATCTGCATTCATTCATGTAAGGCAAAATGGCTGGGTGAGGAGT	575400	59
SEQ ID NO:159	1148T	23-428	DM FAM FRET - WT	FAM-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:160	1148T	23-755	ER24 Red FRET w/FAM stem - MT	Z35-tct-X-agg-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-ggc-ct-hex		
SEQ ID NO:161	711+1G>T	1453-08-06	Invader (+2S)	GCTTTCAGATGTATATTTTAAACAATAGTGCCTAAAGATTAATCAATAGGTACAT	615700	61
SEQ ID NO:162	711+1G>T	1791-08-03	WT DM	CGCGCGAGGACTTCATCAAAATTTGTTTCAG-hex	282800	34

FIGURE 3

SEQ ID NO	Mutation	Oligo Name	Oligo type	Sequence (5'-3')	ε260 M-1 cm-1	Length
SEQ ID NO:163	711+1G>T	1791-06-02	Mt ER38	AGGCCACGGACGAATTCATCAAAATGTTTCAGG-hex	325500	37
SEQ ID NO:164	711+1G>T	1614-56-03	WT ST	AAGCCCTGAACAAATTTGATGAATATGTACCTATTGATTAATCTTTTAGGCACTATTGTTATAAATATACAACCTGGAAG	882700	88
SEQ ID NO:165	711+1G>T	1614-56-07	Mt ST	AACAACCTGAACAAATTTGATGAATATGTACCTATTGATTAATCTTTTAGGCACTATTGTTATAAATATACAACCTGGAAG	880100	88
SEQ ID NO:166	711+1G>T	23-204	ER38 FAM FRET -MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-dt-hex		
SEQ ID NO:167	711+1G>T	23-746	DM Red FRET w/FAMstem-WT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-dt-hex		
SEQ ID NO:168	S549N	1614-62-01	Invader	TCACCTTGCTAAAGAAATTCCTTGCTGCTGACCTCCAA	352200	38
SEQ ID NO:169	S549N	1614-62-04	WT probe (ER24)	ACGGACGCGGAGCTCAGTGTGATTCCACC-hex	275200	33
SEQ ID NO:170	S549N	1614-62-06	MT probe (DM)	CGGCCGAGGTTCACTGTGATTCCACC-hex	247800	31
SEQ ID NO:171	S549N	1614-62-03	WT ST	AGAAAGTGGAATCACACTGAGTGAGGTCAGAGCAACGACGACAAAGAAATTTCTTTAGCAAGGTTGAAT	622500	61
SEQ ID NO:172	S549N	1614-62-07	MT ST	AGAAAGTGGAATCACACTGAGTGAGGTCAGAGCAACGACGACAAAGAAATTTCTTTAGCAAGGTTGAAT	623800	61
SEQ ID NO:173	S549N	23-755	ER24 Red FRET w/FAM stem - WT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:174	S549N	23-428	DM FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-ggc-ggc-hex		
SEQ ID NO:175	D1152H	1645-17-10	Invader	TGAGTACATTGCATGGCTGTAAACTCCAGCATAT	351800	36
SEQ ID NO:176	D1152H	1733-73-02	WT probe (ER24)	ACGGACGCGGAGGATGTGGTAGCTTGGA-hex	303500	34
SEQ ID NO:177	D1152H	1733-73-03	Mt probe (ER38)	AGCCACGGGACGCATGTGGATAGCTTGGA-hex	296300	34
SEQ ID NO:178	D1152H	1645-17-12	WT ST	AGACTTACCAAGTATCCACATGTATGCTGGAGTTACAGCCCACTGCAATGTACTCATGT	581000	61
SEQ ID NO:179	D1152H	1645-17-16	MT ST	AGACTTACCAAGTATCCACATGTATGCTGGAGTTACAGCCCACTGCAATGTACTCATGT	584500	61
SEQ ID NO:180	D1152H	23-210	ER24 FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-gt-hex		
SEQ ID NO:181	D1152H	23-752	ER38 Red FRET w/FAM stem - MT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-ggt-ggc-dt-hex		
SEQ ID NO:182	3905insT	1614-72-12	Invader	GATCTGGATTCTCTCTCAGTGTTCAGTAGTCTCAT	334000	
SEQ ID NO:183	3905insT	1687-91-02	WT probe (ER38)	AGCCACGGACGCAAAAAGCTGATAACAAGTACT	352100	
SEQ ID NO:184	3905insT	1614-72-07	MT probe (DM)	CGCCGCGAGGAAAGAGCTGATAACAAGTACT	337300	
SEQ ID NO:185	3905insT	1645-07-01	WT ST	CAGGGAAGAGTACTTTGTTATCAGCTTTTTTGAGACTACTGAACACTGAAGGAGAAATCCAGATCGATGG	698200	
SEQ ID NO:186	3905insT	1645-07-02	MT ST	CAGGGAAGAGTACTTTGTTATCAGCTTTTTTGAGACTACTGAACACTGAAGGAGAAATCCAGATCGATGG	706300	
SEQ ID NO:187	3905insT	23-428	DM FAM FRET -MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-ggc-ggc-hex		
SEQ ID NO:188	3905insT	23-752	ER38 Red w/FAM stem-WT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-ggc-ggc-hex		
SEQ ID NO:189	Y1092X C>G	1645-14-10	Invader	CCACAAAGCTCTGTAATTTACATACTGCCAACTGGTCTTGTA	404700	
SEQ ID NO:190	Y1092X C>G	1733-75-01	WT probe DM	CGCCGCGAGGCGTCAACACTGC	214800	
SEQ ID NO:191	Y1092X C>G	1645-14-17	MT probe ER24	ACGGACGCGGAGGCTGTCAACACTGCG	258400	
SEQ ID NO:192	Y1092X C>G	1645-14-12	WT ST	CCAGCCGAGTGTGACAGCTACAAGAACCCAGTTGGCAGTATGTAATTCAGAGCTTTGTGGAAT	638100	
SEQ ID NO:193	Y1092X C>G	1645-14-16	MT ST	CCAGCCGAGTGTGACAGCTACAAGAACCCAGTTGGCAGTATGTAATTCAGAGCTTTGTGGAAT	631400	
SEQ ID NO:194	Y1092X C>G	23-428	DM FAM FRET-WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-ggc-ggc-hex		
SEQ ID NO:195	Y1092X C>G	23-394	ER24 Red FRET-MT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-agg-aga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:196	3849+4A>G	1645-13-10	Invader AS	ACATTCTCTTCAATAAGTCTCTGGCCAGAGGGTGT	340400	
SEQ ID NO:197	3849+4A>G	1645-13-14	WT ER38	AGCCACGGAGGAGATTTGAACACTGCTTG	291200	
SEQ ID NO:198	3849+4A>G	1645-13-15	Mt DM	CGCCCGGAGGAGATTTGAACACTGCTTG	261300	
SEQ ID NO:199	3849+4A>G	1645-13-12	WT ST	AAAGCAAGCAGTGTCAAAATCTCAACCTCTGGCCAGGACATTGAGAAGGAAATGTTCT	587100	
SEQ ID NO:200	3849+4A>G	1645-13-16	Mt ST	AAAGCAAGCAGTGTCAAAATCTCAACCTCTGGCCAGGACATTGAGAAGGAAATGTTCT	586200	
SEQ ID NO:201	3849+4A>G	23-205	DM Red FRET	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-dt-hex		
SEQ ID NO:202	3849+4A>G	23-204	ER38 FAM FRET	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-dt-hex		
SEQ ID NO:203	3878delA	1614-61-01	Invader S	CCTTCAGTGTTCAGTAGTCTCAAAAAAGCTGATAACAAGTACTCTTCT	483100	
SEQ ID NO:204	3878delA	1614-61-02	WT DM	CGCCCGGAGGCTCAGTCCAGTCTTCTCCCV	237300	
SEQ ID NO:205	3878delA	1614-61-09	Mt ER38	AGCCACGGAGCGGATCGACTGACTTCTCCCV	256600	
SEQ ID NO:206	3878delA	1614-61-03	WT ST	CTCTTGGGAAGAACTGGATCAGGGAAGAGTACTTTGATCAGCTTTTTTGAGACTACTGAACACTGAAGGAG	719600	
SEQ ID NO:207	3878delA	1614-61-07	MT ST	CTCTTGGGAAGAACTGGATCAGGGAAGAGTACTTTGATCAGCTTTTTTGAGACTACTGAACACTGAAGGAG	706800	
SEQ ID NO:208	3878delA	23-428	DM FAM FRET WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-ttc-ggc-ggc-hex		
SEQ ID NO:209	3878delA	23-752	ER38 Red w/FAM stem MT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-dt-hex		
SEQ ID NO:210	3878delA	1645-12-01	Invader	TGGTGCAGGCATTAATCCAGAAAATCT	275700	
SEQ ID NO:211	Q493X	1767-83-03	WT ER24	ACGGACGCGGAGGAGAACAGAAATGAAATCT	322300	
SEQ ID NO:212	Q493X	1645-12-09	WT ER38	AGCCACGGACGGAAGCAACAGAAATGAAATCTCC	343700	
SEQ ID NO:213	Q493X	1645-12-03	WT ST	CAGTGGGAAGAAATTCATCTGTTCTCAGTTTCTCGGATATGCTCGGCAACCATTT	518600	
SEQ ID NO:214	Q493X	1645-12-07	MT ST	CAGTGGGAAGAAATTCATCTGTTCTTAGTTTCTCGGATATGCTCGGCAACCATTT	520100	
SEQ ID NO:215	Q493X	23-210	ER24 FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-tcc-ggc-tcc-gt-hex		
SEQ ID NO:216	Q493X	23-752	ER38 Red FRET w/FAM stem - MT	235-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-dt-hex		

FIGURE 3

SEQ ID NO	Mutation	Oligo Name	Oligo type	Sequence (5'-3')	ε260 M-1 cm-1	Length
SEQ ID NO:217	G551D	1614-52-10	Invasor AS	TGGAGAAAGGTGGAATCACACTGAGTGAGT	308600	
SEQ ID NO:218	G551D	1711-01-01	WT DM	CGCCGAGGAGGTCAACGAGCAAGAAATTTV	274400	
SEQ ID NO:219	G551D	1711-01-06	MT ER38	AGCCACGAGGATCAACGAGCAAGAAATTTV	311200	
SEQ ID NO:220	G551D	1614-52-12	WT ST	CTAAGAAATCTTGCTGTTGACCTGACCTCAGTGTGATCCACCTTCTCCAAAGT	514900	
SEQ ID NO:221	G551D	1614-52-16	MT ST	CTAAGAAATCTTGCTGTTGATCTCAGTGTGATCCACCTTCTCCAAAGT	516800	
SEQ ID NO:222	G551D	23-746	DM Red FRET w/FAM stem WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:223	G551D	23-204	ER38 FAM FRET	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:224	R553X	1645-72-01	Invasor +2S	GCTAGACCAATAATGATTCACCTTGCTAAAGAAATCTTGCTG	456700	
SEQ ID NO:225	R553X	1667-08-02	WT ER38	AGCCACGAGGAGGCTGACCTCCACTCA	260300	
SEQ ID NO:226	R553X	1453-07-01	MT DM	CGCCGAGGAGGATGACCTCCACTCAGT	263700	
SEQ ID NO:227	R553X	1614-53-03	WT ST	TCACACTGAGTGAGGTCAACGAGCAAGAAATTTCTTTAGCAAGGTGAATACTAATTTATGGTCTAGCAAGCT	727300	
SEQ ID NO:228	R553X	1614-53-07	MT ST	TCACACTGAGTGAGGTCAATGAGCAAGAAATTTCTTTAGCAAGGTGAATACTAATTTATGGTCTAGCAAGCT	728600	
SEQ ID NO:229	R553X	23-428	DM FAM FRET MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:230	R553X	23-752	ER38 Red FRET w/FAM stem WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:231	R1162X	1618-52-01	Invasor	GTTTACCTTCTTGTCGATGTCAATGAACTTAAAGACTCT	424200	
SEQ ID NO:232	R1162X	1618-52-03	WT probe (ER24)	AGGACGCGGAGGCTCACAGATCGCV	284000	
SEQ ID NO:233	R1162X	1618-52-08	MT probe (ER38)	AGCCACGAGGAGGCTCACAGATCGCV	283800	
SEQ ID NO:234	R1162X	1618-52-05	WT ST	AGATCGATCTGTGAGCGAGTCTTTAAGTTTATGATGACATGCCAACAGAGGTAAGCT	574600	
SEQ ID NO:235	R1162X	1618-52-09	MT ST	AGATCGATCTGTGAGCTGAGTCTTTAAGTTTATGATGACATGCCAACAGAGGTAAGCT	574900	
SEQ ID NO:236	R1162X	23-204	ER38 FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:237	R1162X	23-755	ER24 Red FRET w/FAM - WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:238	S549R (A>C)	1614-63-01	Invasor	CACCTTGTAAAGAAATTTCTTGCTGTGACCTCCACC	344400	
SEQ ID NO:239	S549R (A>C)	1614-63-02	WT probe (DM)	CGCCGAGGAGGCTGATCCACTV	247500	
SEQ ID NO:240	S549R (A>C)	1733-78-02	MT probe (ER38)	AGCCACGAGGAGGCTGATCCACTV	259500	
SEQ ID NO:241	S549R (A>C)	1614-63-03	WT ST	GAGAAGGTGGAATCACACTGAGTGGAGGTCAACGAGCAAGAAATTTCTTTAGCAAGGTGAAT	624200	
SEQ ID NO:242	S549R (A>C)	1614-63-07	MT ST	GAGAAGGTGGAATCACACTGCGTGGAGGTCAACGAGCAAGAAATTTCTTTAGCAAGGTGAAT	617600	
SEQ ID NO:243	S549R (A>C)	23-428	DM FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:244	S549R (A>C)	23-211	ER38 Red FRET - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:245	F508C	1618-70-10	Invasor AS	CAGTTTCTGATGTGCTGCGACCATTAAGAAATATCATCTC	450800	
SEQ ID NO:246	F508C	1791-37-01	WT DM	CGCCGAGGAGGCTGATCCACTV	291900	
SEQ ID NO:247	F508C	1828-07-05	MT ER38	AGCCACGAGGAGGCTGATCCACTV	287800	
SEQ ID NO:248	F508C	1618-70-14	WT ST	TATATCATATAGAAACACCAAGATGATATTTCTTTAAGTGGCAGGCAATCCAGGAAACTGAGT	736000	
SEQ ID NO:249	F508C	1618-70-18	MT ST	TATATCATATAGAAACACCAAGATGATATTTCTTTAAGTGGCAGGCAATCCAGGAAACTGAGT	731600	
SEQ ID NO:250	F508C	23-428	DM FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:251	F508C	23-211	ER38 Red FRET Mt	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:252	Y1092X C>A	1645-14-10	Invasor	CCACAAAGCTCTGAAATTTACATAGCTGCAACTGCTTGTAT	404700	
SEQ ID NO:253	Y1092X C>A	1645-14-11	WT probe -DM	CGCCGAGGAGGCTGCAACACTGCG	225400	
SEQ ID NO:254	Y1092X C>A	1645-15-08	MT probe -ER38	AGCCACGAGGAGGCTGCAACACTGCG	258400	
SEQ ID NO:255	Y1092X C>A	1645-14-12	WT ST	CCAGCGAGTGTGACAGGTACAAGCAACCAAGTGGCAGTGTAAATTCAGAGCTTTGTGGAAT	636100	
SEQ ID NO:256	Y1092X C>A	1645-15-06	MT ST	CCAGCGAGTGTGACAGTACAAGCAACCAAGTGGCAGTGTAAATTCAGAGCTTTGTGGAAT	634100	
SEQ ID NO:257	Y1092X C>A	23-428	DM FAM FRET - WT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:258	Y1092X C>A	23-752	ER38 Red FRET w/FAM - MT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:259	S549R (T>G)	1614-64-01	Invasor	AGTTATACCTTGTGTAAGAAATTTCTTGCTGTTGACCTCCT	397300	
SEQ ID NO:260	S549R (T>G)	1778-64-06	WT ER38	AGCCACGAGGAGGCTGATCCACTV	285300	
SEQ ID NO:261	S549R (T>G)	1614-64-06	MT probe (DM)	CGCCGAGGAGGCTGATCCACTV	245100	
SEQ ID NO:262	S549R (T>G)	1614-64-03	WT ST	GAAAGTGGAAATCACACTGAGGAGGAGGTCACAGGCAAGAAATTTCTTTAGCAAGGTGAATACTAAT	675300	
SEQ ID NO:263	S549R (T>G)	1614-64-07	MT ST	GAAAGTGGAAATCACACTGAGGAGGAGGTCACAGGCAAGAAATTTCTTTAGCAAGGTGAATACTAAT	676700	
SEQ ID NO:264	S549R (T>G)	23-211	ER38 Red FRET - WT	Z35-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:265	S549R (T>G)	23-428	DM FAM FRET - MT	FAM-tdt-X-agg-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-gcg-hex		
SEQ ID NO:266	del1507	1614-58-01	Invasor S	CATGCTTGTGACGCTTCTGATCTATATTCATATGATAGGAAACACCAAT	493900	
SEQ ID NO:267	del1507	1645-58-01	WT ER24	ACGGACGCGGAGGATGATATTTCTTTAATGGTG	339900	
SEQ ID NO:268	del1507	1645-58-04	MT ER38	AGCCACGAGGAGGATATTTCTTTAATGGTG	314500	
SEQ ID NO:269	del1507	1614-58-03	WT ST	ATGCTGGCACCATTAAAGAAATATCATCTTTGTTGTTTCTTATGATGAATATAGACAGGCTCATCAAGCATGCC	814200	
SEQ ID NO:270	del1507	1614-58-07	MT ST	ATGCTGGCACCATTAAAGAAATATCATCTTTGTTGTTTCTTATGATGAATATAGACAGGCTCATCAAGCATGCC	785500	

FIGURE 3

SEQ ID NO	Mutation	Oligo Name	Oligo type	Sequence (5'-3')	ε260 M-1 cm-1	Length
SEQ ID NO:271	del1507	23-204	ER38 FAM FRET MT	FAM-tct-X-agc-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:272	del1507	23-755	ER24 Redw/FAM FRET WT	Z35-tct-X-agc-cgg-ttt-tcc-ggc-tga-gac-tcc-gcg-tcc-gt-hex		
SEQ ID NO:273	5T Variant (5T)	16924810	Probe (DM)	CGGCCGAGGGTTTTTAACAGGGGATTTGGG	286900	
SEQ ID NO:274	5T Variant (7T)	16925804	Probe (DM)	CGGCCGAGGGTTTTTAACAGGGGATTTGG	293000	
SEQ ID NO:275	5T Variant (9T)	16929202	Probe (DM)	CGGCCGAGGGTTTTTTTTTACCAGGGATTGGGA	338300	
SEQ ID NO:276	5T Variant	16027813	Invader	CTCATCTTTATTTGATGTGTGTGTGTGTGTGA	353300	
SEQ ID NO:277	5T Variant	16027814	Invader	CTCATCTTTATTTGATGTGTGTGTGTGTGTGTGA	372100	
SEQ ID NO:278	5T Variant	16027815	Invader	CTCATCTTTATTTGATGTGTGTGTGTGTGTGTGTGA	390900	
SEQ ID NO:279	5T Variant	16027816	Invader	CTCATCTTTATTTGATGTGTGTGTGTGTGTGTGTGTGA	409700	
SEQ ID NO:280	5T Variant	16060401	Invader	CTCATCTTTATTTGATGTGTGTGTGTGTGTGTGA	334500	
SEQ ID NO:281	5T Variant	23-428	DM FAM FRET	FAM-tct-X-agc-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gcg-hex		
SEQ ID NO:282	Y122X	1809-46-01	Invader	GTGTCTTCACAATAAAGAGAGGCATAGCCTATGCCTAA	403500	
SEQ ID NO:283	Y122X	1827-34-03	WT ER24	ACGGACGGGAGGATAAATCGCGGATAGAG	309000	
SEQ ID NO:284	Y122X	1827-34-12	MT ER38	AGGCCAGCGAGCGTTAAATCGCGGATAGAG	295100	
SEQ ID NO:285	Y122X	1645-16-03	WT ST	GGAACGCTCTATCGCGATTATCTAGGCATAGGCTTATGCTCTCTTTATGTGAGGACACTGT	609400	
SEQ ID NO:286	Y122X	1645-16-07	WT ST	GGAACGCTCTATCGCGATTAACTAGGCATAGGCTTATGCTCTCTTTATGTGAGGACACTGT	612300	
SEQ ID NO:287	Y122X	23-394	ER24 Red FRET -WT	Z35-tct-X-tgg-ggc-ttt-tgg-cgg-aga-gac-tcc-gcg-tcc-hex		
SEQ ID NO:288	Y122X	23-204	ER38 FAM FRET - MT	FAM-tct-X-agc-cgg-ttt-tcc-ggc-tga-gac-gtc-cgt-ggc-ct-hex		
SEQ ID NO:289	1898+1 G>A	1353-02-02	Invader	GACTCTCTCTTTTGATACCTTAGATGTTTAAACAGAAAAGAAATATTTGAAAGT		
SEQ ID NO:290	1898+1 G>A	1801-65-03	WT ER24	ACGGACGGCGAGTATGTCTTTGAAATACCTTACT		
SEQ ID NO:291	1898+1 G>A	1801-65-08	MT DM	CGGCCGAGGATATGTCTTTGAAATACCTTACTT		
SEQ ID NO:292	1898+1 G>A	1614-59-12	WT ST	AGCATATAAGTAAAGTATTCAAAGAACATACCTTTCAAAATATCTTTCTGTTAAACATCTAGGTATCCAAAAGGAGT		
SEQ ID NO:293	1898+1 G>A	1614-59-16	MT ST	AGCATATAAGTAAAGTATTCAAAGAACATATCTTTCAAAATATCTTTCTGTTAAACATCTAGGTATCCAAAAGGAGT		
SEQ ID NO:294	1898+1 G>A	23-210	ER24 FAM FRET - WT	FAM-tct-X-agc-cgg-ttt-tcc-ggc-tga-gac-tcc-gcg-tcc-gt-hex		
SEQ ID NO:295	1898+1 G>A	23-746	DM Red FRET w/FAM- MT	Z35-tct-X-agc-cgg-ttt-tcc-ggc-tga-gac-ctc-ggc-gcg-hex		

X = Quencher = Z28

9= 2' O-methyl modified base

FIGURE 4

SEQ ID NO	Allele	Part No	1614-48-17	Mt Probe (ER24)	ACGGACGCGGAGACTCTATCGCGATTATCTAGV
SEQ ID NO:417	R117H				
SEQ ID NO:418	R117H	1614-48-10		invader	AGAAATCATAGCTTCCTATGACCCCGGATAACAAGGAGAACT
SEQ ID NO:419	W1282X	1618-42-04		Wt probe (ER38)	AGGCCACGGACGCCACTGTTGCAAAAGTTATT
SEQ ID NO:420	711+1G>T	1791-06-04		Wt ER24	ACGGACGCGGAGACTTCATCAAAATTTGTTTCAG
SEQ ID NO:421	711+1G>T	1510-82-03		Mt ER24	ACGGACGCGGAGAAATTCATCAAAATTTGTTTCAGG
SEQ ID NO:422	711+1G>T	1614-56-05		Wt probe (ER38)	AGGCCACGGACGCTTCATCAAAATTTGTTTCAGGTV
SEQ ID NO:423	711+1G>T	1614-56-06		Mt probe (DM)	CGGCCGAGGATTCATCAAAATTTGTTTCAGGTV
SEQ ID NO:424	711+1G>T	1614-56-08		Mt Probe (ER24)	ACGGACGCGGAGATTCATCAAAATTTGTTTCAGGTV
SEQ ID NO:425	711+1G>T	1614-56-01		invader	CCTTCCAGTTGTATAATTTATAACAATAGTGCCTAAAGATT AAATCAATAGGTACATAT
SEQ ID NO:426	394delTT			invader	CCCTTCGGCGATGTTTTTCTGGAGATTTATGTTCTATGGAA
SEQ ID NO:427	394delTT			Tt	
SEQ ID NO:428	394delTT			WT probe	AGGCCACGGACGCTTTTATATTTAGGGTAAGGATCTHex
				MT probe	AAGCACGCAGCACCTTTATATTTAGGGTAAGGATCTCHex
SEQ ID NO:429	del507	1353-01-01		Invader S	GCTTTGATGACGCTTCTGTATCTATATTCATCATAGGAAACAC
SEQ ID NO:430	del507	1556-23-09		Wt DM	CAAT
SEQ ID NO:431	del507	1826-61-01		Wt DM	CGGCCGAGGAGATGATATTTCTTTAATGGTG
SEQ ID NO:432	del507	1556-23-10		WT ER24	CGGCCGAGGAGATGATATTTCTTTAATGGT
SEQ ID NO:433	del507	1826-61-02		WT ER24	ACGGACGCGGAGAGATGATATTTCTTTAATGGTG
SEQ ID NO:434	del507	1556-23-11		WT ER38	ACGGACGCGGAGAGATGATATTTCTTTAATGGT
SEQ ID NO:435	del507	1826-61-03		WT ER38	AGGCCACGGACGAGATGATATTTCTTTAATGGTG
SEQ ID NO:436	del507	1353-01-03		Mt DM	AGGCCACGGACGAGATGATATTTCTTTAATGGT
SEQ ID NO:437	del507	1826-61-04		Mt DM	CGGCCGAGGAGATATTTCTTTAATGGTCCC
SEQ ID NO:438	del507	1353-55-01		Mt ER24	CGGCCGAGGAGATATTTCTTTAATGGTGC
SEQ ID NO:439	del507	1826-61-05		Mt ER24	ACGGACGCGGAGAGATATTTCTTTAATGGTGCCA
SEQ ID NO:440	del507	1826-61-06		Mt ER24	ACGGACGCGGAGAGATATTTCTTTAATGGTGCC
SEQ ID NO:441	del507	1353-55-02		Mt ER38	ACGGACGCGGAGAGATATTTCTTTAATGGTGC
SEQ ID NO:442	del507	1826-61-07		Mt ER38	AGGCCACGGACGAGATATTTCTTTAATGGTGCCA
SEQ ID NO:443	del507	1826-61-08		Mt ER38	AGGCCACGGACGAGATATTTCTTTAATGGTGCC
SEQ ID NO:444	5T Variant (5T)	18312001		Probe (ER24)	AGGCCACGGACGAGATATTTCTTTAATGGTGC
SEQ ID NO:445	5T Variant (5T)	18312002		Probe (ER38)	ACGGACGCGGAGGTTTTTAACAGGGATTTGGG
SEQ ID NO:446	5T Variant (5T)	16027817		Probe (DM)	AGGCCACGGACGCGTTTTTAACAGGGATTTGGG
SEQ ID NO:447	5T Variant (5T)	16027818		Probe (ER24)	CGGCCGAGGGTTTTTAACAGGGATTTGGG
SEQ ID NO:448	5T Variant (5T)	16027819		Probe (ER38)	ACGGACGCGGAGGTTTTTAACAGGGATTTGGG
					AGGCCACGGACGCGTTTTTAACAGGGATTTGGG

FIGURE 4

SEQ ID NO:449	5T Variant (7T)	18312003	Probe (ER24)	ACGACGCGGAGGTTTTTTTAAACAGGATTG
SEQ ID NO:450	5T Variant (7T)	18312004	Probe (ER38)	AGCCACGACGCGTTTTTTTAAACAGGATTG
SEQ ID NO:451	5T Variant (7T)	16924813	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG
SEQ ID NO:452	5T Variant (7T)	16924814	Probe (ER24)	ACGACGCGGAGGTTTTTTTAAACAGGATTG
SEQ ID NO:453	5T Variant (7T)	16924815	Probe (ER38)	AGCCACGACGCGTTTTTTTAAACAGGATTG
SEQ ID NO:454	5T Variant (7T)	16027820	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTG
SEQ ID NO:455	5T Variant (7T)	16027821	Probe (ER24)	ACGACGCGGAGGTTTTTTTAAACAGGATTG
SEQ ID NO:456	5T Variant (7T)	16027822	Probe (ER38)	AGCCACGACGCGTTTTTTTAAACAGGATTG
SEQ ID NO:457	5T Variant (9T)	18312005	Probe (ER24)	ACGACGCGGAGG9TTTTTTTACCAGGATTGGGA*
SEQ ID NO:458	5T Variant (9T)	18312006	Probe (ER38)	AGCCACGACGCG9TTTTTTTACCAGGATTGGGA*
SEQ ID NO:459	5T Variant (9T)	16929204	Probe (DM)	CGGCCGAGGG9TTTTTTTAAACAGGAAATTGGGA*
SEQ ID NO:460	5T Variant (9T)	18312007	Probe (ER24)	ACGACGCGGAGG9TTTTTTTAAACAGGAAATTGGGA*
SEQ ID NO:461	5T Variant (9T)	18312008	Probe (ER38)	AGCCACGACGCG9TTTTTTTAAACAGGAAATTGGGA*
SEQ ID NO:462	5T Variant (9T)	17337201	Probe (DM)	CGGCCGAGGG9TTTTTTTACCAGGATTGGGA*
SEQ ID NO:463	5T Variant (9T)	18312009	Probe (ER24)	ACGACGCGGAGG9TTTTTTTACCAGGATTGGGA*
SEQ ID NO:464	5T Variant (9T)	18312010	Probe (ER38)	AGCCACGACGCG9TTTTTTTACCAGGATTGGGA*
SEQ ID NO:465	5T Variant (9T)	17337202	Probe (DM)	CGGCCGAGGG9TTTTTTTAAACAGGAAATTGGGA*
SEQ ID NO:466	5T Variant (9T)	18312011	Probe (ER24)	ACGACGCGGAGG9TTTTTTTAAACAGGAAATTGGGA*
SEQ ID NO:467	5T Variant (9T)	18312012	Probe (ER38)	AGCCACGACGCG9TTTTTTTAAACAGGAAATTGGGA*
SEQ ID NO:468	5T Variant (9T)	16027823	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG
SEQ ID NO:469	5T Variant (9T)	16027824	Probe (ER24)	ACGACGCGGAGGTTTTTTTAAACAGGATTGGG
SEQ ID NO:470	5T Variant (9T)	16027825	Probe (ER38)	AGCCACGACGCGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:471	5T Variant (9T)	16925101	Probe (DM)	CGGCCGAGGGTTATTTTAAACAGGATTGGG*
SEQ ID NO:472	5T Variant (9T)	16925102	Probe (DM)	CGGCCGAGGGTTATTTTAAACAGGATTGGG*
SEQ ID NO:473	5T Variant (9T)	16925103	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:474	5T Variant (9T)	16925104	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:475	5T Variant (9T)	16925105	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:476	5T Variant (9T)	16925106	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:477	5T Variant (9T)	16925107	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:478	5T Variant (9T)	16925108	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:479	5T Variant (9T)	16925109	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:480	5T Variant (9T)	16925110	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:481	5T Variant (9T)	16925111	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:482	5T Variant (9T)	16925112	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:483	5T Variant (9T)	16925113	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:484	5T Variant (9T)	16925114	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*
SEQ ID NO:485	5T Variant (9T)	16925115	Probe (DM)	CGGCCGAGGGTTTTTTTAAACAGGATTGGG*

FIGURE 4

SEQ ID NO:486	5T Variant (9T)	16925116	Probe (DM)	CGCGCCGAGGGGTTTTTTTTTAAACAGGGGATATGGGG*
SEQ ID NO:487	5T Variant (9T)	16925117	Probe (DM)	CGCGCCGAGGGGTTTTTTTTTAAACAGGGGATTAGGGG*
SEQ ID NO:488	5T Variant (9T)	16925118	Probe (DM)	CGCGCCGAGGGGTTTTTTTTTAAACAGGGGATTTAGGG*
SEQ ID NO:489	5T Variant (9T)	16925119	Probe (DM)	CGCGCCGAGGGGTTTTTTTTTAAACAGGGGATTTGAGG*
SEQ ID NO:490	5T Variant (9T)	16925120	Probe (DM)	CGCGCCGAGGGGTTTTTTTTTAAACAGGGGATTTGGAG*
SEQ ID NO:491	5T Variant (9T)	16927101	Probe (DM)	CGCGCCGAGGGG9TTTTTTTTTAAACAGGGGATTTGGGG*
SEQ ID NO:492	5T Variant (9T)	16927102	Probe (DM)	CGCGCCGAGGGG9TTTTTTTTTAAACAGGGGATTTGGGG*
SEQ ID NO:493	5T Variant (9T)	16929201	Probe (DM)	CGCGCCGAGGGG9TTTTTTTTTACCAGGGGATTTGGGG*
SEQ ID NO:494	5T Variant (9T)	16929203	Probe (DM)	CGCGCCGAGGGG9TTTTTTTTTAAACAGGGGAATTGGGG*

9 = 2' O-methyl modified base

FIGURE 5

SEQ ID NO		forward primer
SEQ ID NO: 394	cftr exon 3	TGGTCCCACCTTTTATTCTTTTGCAGA
SEQ ID NO: 395	cftr exon 4	AAGTCACCAAAGCAGTACAGCC
SEQ ID NO: 396	cftr exon 5	GCTGTCAAGCCGTGTTCTAGATAAA
SEQ ID NO: 397	cftr exon 7	CGGAAGGCAGCCTATGTGAGA
SEQ ID NO: 398	cftr exon 9	CATGGGCCATGTGCTTTTCAAAC
SEQ ID NO: 399	cftr exon 9-1	CATGGGCCATGTGCTTTTCAAAC
SEQ ID NO: 400	cftr exon 9-2	CTTCTTGGTACTCCTGTCCTGAAAGA
SEQ ID NO: 401	cftr exon 10	ATTATGGGAGAACTGGAGCCTTCA
SEQ ID NO: 402	cftr exon 11	GATTACATTAGAAGGAAGATGTGCCTTTCAA
SEQ ID NO: 403	cftr exon 12	TAAGGCCAAATCATCTACACTAGATGACCA
SEQ ID NO: 404	cftr exon 13	TAACCTGAGACCTTACACCGTTTCTCA
SEQ ID NO: 405	cftr exon 14B	ATGGGAGGAATAGGTGAAGATGTTAGAA
SEQ ID NO: 406	cftr exon 16	TCTGAATGCGTCTACTGTGATCCA
SEQ ID NO: 407	cftr exon 17A	CCTGCACAATGTGCACATGTACC
SEQ ID NO: 408	cftr exon 17B	GGACTATGGACACTTCGTGCC
SEQ ID NO: 409	cftr exon 18	GGAGAAGGAAGAGTTGGTATTATCCTGAC
SEQ ID NO: 410	cftr exon 19	GCATCAAACCTAATTGTGAAATTGTCTGCC
SEQ ID NO: 411	cftr exon 19-1	GCATCAAACCTAATTGTGAAATTGTCTGCC
SEQ ID NO: 412	cftr exon 19-2	GAAGGTGGAAATGCCATATTAGAGAACA
SEQ ID NO: 413	cftr exon 20	GTACCTATATGTCACAGAAGTGATCCCA
SEQ ID NO: 414	cftr exon 21	GATTAGAAAAATGTTCAAGGGACTCCA
SEQ ID NO: 415	cftr 3849+10kb	CAGTTGACTTGTCTCTTGATTTCTGGA
SEQ ID NO: 416	cftr exon 17A-2	CCTCGACAATGTGCACATGTACC

		reverse primer
SEQ ID NO: 495	cftr exon 3	ACCTATTCACCAGATTTTCGTAGTCTTTTCA
SEQ ID NO: 496	cftr exon 4	TGTACCAGCTCACTACCTAATTTATGACA
SEQ ID NO: 497	cftr exon 5	GAGCTGAGCAAGACTTAACCACTAATTAC
SEQ ID NO: 498	cftr exon 7	GTGAACATTCTAGTATTAGCTGGCAAC
SEQ ID NO: 499	cftr exon 9	CTCCAAAAATACCTTCCAGCACTACAAA
SEQ ID NO: 500	cftr exon 9-1	GAAATTACTGAAGAAGAGGCTGTCATCAC
SEQ ID NO: 501	cftr exon 9-2	CTCCAAAAATACCTTCCAGCACTACAAA
SEQ ID NO: 502	cftr exon 10	GACTAACCGATTGAATATGGAGCCAAA
SEQ ID NO: 503	cftr exon 11	CTTAAATGTGATTCTTAACCCACTAGCCA
SEQ ID NO: 504	cftr exon 12	GAGGTAAAATGCAATCTATGATGGGACA
SEQ ID NO: 505	cftr exon 13	TAAGGGAGTCTTTTGCACAATGGAAAA
SEQ ID NO: 506	cftr exon 14B	ACCTCACCCAACTAATGGTCATCA
SEQ ID NO: 507	cftr exon 16	TAGACAGGACTTCAACCCCTCAATCA
SEQ ID NO: 508	cftr exon 17A	GAGTATCGCACATTCACTGTCATACC
SEQ ID NO: 509	cftr exon 17B	AAGGTAACAGCAATGAAGAAGATGACAAA
SEQ ID NO: 510	cftr exon 18	TAATGACAGATACACAGTGACCCTCAA
SEQ ID NO: 511	cftr exon 19	GCTTCAGGCTACTGGGATTAC
SEQ ID NO: 512	cftr exon 19-1	GTCATCTTTCTTCACGTGTGAATTCTCAA
SEQ ID NO: 513	cftr exon 19-2	GCTTCAGGCTACTGGGATTAC
SEQ ID NO: 514	cftr exon 20	TTCTGGCTAAGTCCTTTTGCTCAC
SEQ ID NO: 515	cftr exon 21	CATTTTCAGTTAGCAGCCTTACCTCA
SEQ ID NO: 516	cftr 3849+10kb	TCCTCCCTGAGAATGTTGGATCAA
SEQ ID NO: 517	Vs1 Int std F	TGATGGTGGTATGTTTTTCAGGCTAGA
SEQ ID NO: 518	Vs1 Int std R	GTTCTCCCCTGTCCAGTTTAAAC

Fig. 6**A**

Mutation	Sample	IC ALLELE	MUT ALLELE	FOZ Ratio
2789+5G>A	26mix	3.94	4.69	1.19
R1162X	29	3.42	2.18	0.62
R347P	15	3.38	4.60	1.36
G85E	21	3.62	2.55	0.70
R560T	9	3.30	2.47	0.75
delI507	1	3.16	1.98	0.63
1898+1G>A	111 A2/8	6.23	2.84	0.46
R117H	30	3.46	1.87	0.54
delF508 homo MT	3	3.44	1.14	0.33
WT gDNA	03-243	3.58	1.06	0.30

Mutation	Sample	IC ALLELE	MUT ALLELE	FOZ Ratio
2184delA plasmid/internal control syn. Target	plasmid/syn. Target	4.67	3.65	0.78

B

Mutation	Sample	IC ALLELE	MUT ALLELE	FOZ Ratio
A455E	8	3.26	2.88	0.88
3659delC	14	3.38	2.36	0.68
N1303K	16	3.92	2.11	0.54
3120+1G>A	6	3.84	2.45	0.64
G551D	20	3.44	2.04	0.59
WT gDNA	03-243	3.74	1.00	0.27

I148T/Internal control	syn. target	4.35	5.08	1.17
1078delT/Internal control	syn. target	4.44	4.97	1.12

C


Mutation	Sample	IC ALLELE	MUT ALLELE	FOZ Ratio
711+1G>T	2	3.95	2.82	0.71
W1282X	19	4.44	2.16	0.49
1717-1G>A	28	4.87	2.19	0.45
3849+10kbC>T	5	3.82	2.48	0.65
WT gDNA	03-243	4.67	1.10	0.24

D

Mutation	Sample	IC ALLELE	MUT ALLELE	FOZ Ratio
621+1G>T	11	4.23	2.05	0.49
G542X	18	3.40	2.83	0.81
R553X	7	4.53	3.27	0.72
R334W	22	3.72	2.79	0.75
WT gDNA	03-243	4.18	1.14	0.27

FIGURE 7

All Sample Extraction Methods




G85E

Legend

- Red = synthetic Mut
- Blue = synthetic WT
- Green = HET
- Black = WT or NTB


All Sample Extraction Methods



A455E

Legend
Red = synthetic Mut
Blue = synthetic WT
Green = HET
Black = WT or NTB

All Sample Extraction Methods



3659 del C

Legend

Red = synthetic Mut.
Blue = synthetic WT
Green = HET
Black = WT or NTB

Figure 8

Analysis of Characterized Samples _ Range Determination




Figure 9A





Figure 9B




Figure 9C




Figure 9D

Genomic C21, 21plex amplified,
tested with individual CF Invader assays or CFVs 1 pools

